

**ACE Stopper size options**

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

**ACE stopper afmetingen**

**medicina** part of HMC Group  
**ACE stopper & dressing**  
 Silicone stopper for antegrade continence enemas (ACE).

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**CAREFULLY READ THE INSTRUCTIONS BEFORE USE.**

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**ACE Stopper**  
**SINGLE PATIENT USE** **STERILE EO**  
 Sterilized using Ethylene Oxide  
 Do not use if packaging is damaged. Keep dry. Keep away from sunlight. Consult instructions for use. Caution. Does not contain Natural Rubber Latex. DEHP. Not made with phthalates, plasticizers or BPA.  
**CAUTION:** FEDERAL LAW (USA) restricts this device to sale by or on the order of a physician. Follow Physicians instructions using this device.  
 Rx only. Caution: This device is for sale by or on order of a physician.

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**SD03 Stopper Dressing**  
**SINGLE PATIENT USE** **STERILE EO**  
 Sterilized using Ethylene Oxide  
 Single Use. Do not re-sterilize. Do not use if packaging is damaged. Keep dry. Keep away from sunlight. Caution. Consult instructions for use.  
**CAUTION:** FEDERAL LAW (USA) restricts this device to sale by or on the order of a physician. Follow Physicians instructions using this device.  
 Rx only. Caution: This device is for sale by or on order of a physician.

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**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Medicina Ltd.**  
 Units 2-4, Rivington View Business Park, Station Road, Bolton BL6 5BN UK  
 CE 1639

**Description** - Medicina ACE stoppers are implant grade silicone stoppers used to seal stoma sites. Soft, comfortable and in a variety of lengths and shaft diameters for different size of the stoma.

- ACE stoppers are indicated for patients with antegrade continence enemas (ACE), sealing the tract or maintaining patency between insertions of the enema catheter.
- ACE stoppers are also indicated for stomas formed for catheter introduction into the bladder such as Mitrofanoff procedures. Note: They can also be used for any stoma of suitable size anywhere in the body.


**⚠ Caution** - The ACE Stopper can be used for up to 30 days of intermittent use, and not used on patients with a history of allergy to silicone. If skin irritation occurs, a small amount of petroleum jelly or a keyhole dressing may form a temporary barrier. Consult Stoma Nurse if irritation persists.

**Placement** - The site should be cleaned and dried before placement. The stopper can be lubricated with water based lubricant if necessary. The stopper is then inserted gently into the tract in the same way as the catheter. When the skin disk is flat to the skin it can be secured with the ACElok dressing (SD03). These dressings are changed every 3-4 days or when they become soiled.

**Removal** - The stopper is removed when the catheter is to be used. The stopper should be cleaned with soap and water, rinsed under the tap and stored in a clean dry place until catheterisation is complete. (Or in line with local policies).

**Maintenance** - Ace stoppers should be replaced after a maximum of 30 days of intermittent use. They are to be used by a single patient only. They should be cleaned with soap and water after each insertion and allowed to dry before storage in a clean dry place. Stoppers should be replaced if there is any damage to the joint between the skin disc and the stopper shaft.

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