

150

300



PEEL OPEN



REF
Catalogue number

CODE

Description
Descrizione

BARCODE GS1-128
GTIN - EXP DATE - LOT

LOT
Batch Code

XXXXXXXX

Date of manufacture

YYYY-MM-GG

Use By

YYYY-MM-GG

Approximate volume. Use for single patient only. Change every 24 hrs.

The device must be used exclusively by healthcare professionals.

Warning: to avoid reflux, place the bag upright and at a level lower than the patient's stomach.

Volume approssimato. Per uso su singolo paziente. Sostituire ogni 24 ore.

Il dispositivo deve essere utilizzato solo da personale sanitario qualificato.

Attenzione: per evitare reflussi, posizionare la sacca in verticale ed a livello più basso rispetto allo stomaco del paziente.

Volume approximatif. Pour une utilisation par patient unique. Changer toutes les 24 heures.

Le dispositif doit être utilisé par personnel de santé qualifié exclusivement.

Attention: pour éviter le reflux, placez le sac à la verticale et à un niveau plus bas que l'estomac du patient.

Ungefähre Füllmenge. Nur für einen Patienten verwenden. Wechseln Sie alle 24 Stunden.

Das Gerät darf nur von medizinischen Fachpersonal verwendet.

Wichtiger Hinweis: Zur Vermeidung von Reflux sollte der Beutel senkrecht und tiefer als der Magen des Patienten platziert werden

El volumen aproximado. Para el uso en un solo paciente. Debe cambiarse cada 24 horas.

El dispositivo tiene que ser utilizado por personal de salud calificado.

Advertencia: Para evitar el reflujo, coloque la bolsa en posición vertical y a un nivel más abajo que el estómago del paciente.

CAUTION: FEDERAL LAW (USA) restricts this device to sale by or on the order of a physician. Follow Physicians instructions using this device.



Do not reuse

STERILE EO

Sterilized using Ethylene Oxide



Do not resterilize

MD

Medical Device

CE

0051

THE DEVICE IS NOT MADE WITH NATURAL RUBBER LATEX



U.S. Federal law restricts this device to sale by or on order of a physician only



Keep Dry



Keep away from sunlight



Do not use if package is damaged

E7309138 REV.02_10/2021



HMC Premedical S.p.A.
Via Bosco, 1/3 - 41037 Mirandola (MO) ITALY
Tel. +39 0535 22704 - Fax. +39 0535 609546

